







Press and hold the tip of
your tongue to the back
two electrodes to move
the cursor up!






Press and hold the tip of
your tongue to the front
two electrodes to move
the cursor down!







Press and hold the tip of
your tongue to the left
electrode to move the
cursor left!







Press and hold the tip of
your tongue to the right
electrode to move the
cursor right!







Flatten your tongue against the array to touch electrodes then remove your tongue from the array to click! Do this twice quickly to double click!



Hold Your tongue against all electrodes for 2 seconds to Right Click!



Hold Your tongue against all electrodes, and without moving your tongue off the array, slide your tongue to a directional electrode to click and drag!